



Genomics Core Lab - Rm N109
Dalhousie University
1348 Summer Street
Halifax, Nova Scotia, B3H 0A8

Contact: Genomics@dal.ca

Sample submission guidelines for 10X Genomics Visium Platform

Prior to sample submission:

Visium runs are time sensitive so please contact Genomics@dal.ca prior to beginning your assay so we can coordinate slide preparation and instrument availability.

Complete the sample submission form.

Provide a cell phone number for a point of contact for any potential issues during assay preparation.

Sample drop off or shipping:

If dropping off samples by hand to the Genomics Core lab, contact Genomics@dal.ca to ensure a seamless handover.

If shipping samples, please only ship on Monday or Tuesday to prevent weekend delays.

For fresh frozen samples place a generous supply of dry ice to ensure dry ice will remain for the duration of the delivery time.

For international clients, we recommend shipping with World Courier for tissues/cells.

Within Canada, we recommend FedEx Next Day Priority services.

Shipping address:

Attn: Mat Nightingale

Genomics Core Lab

1348 Summer Street

Life Science Research Institute

Dalhousie University

Halifax, Nova Scotia, B3H 0A8

Receiving Conditions

Fresh Frozen Samples

RNA quality

- RNA quality assessment should be carried out before placing tissue sections on slides.
- RNA Integrity Number (RIN) should be >4.
- If you'd like us to perform this QC step please supply two 20 µm FFPE sections prior to preparing your slides.

Section placement on slides:

- Place the tissue sections within the allowable area. Place one tissue section per slide unless the sections are small enough to be captured together.
- Ensure that tissue sections and associated OCT compound **do not** overlap with other tissue sections and associated OCT compound if placing multiple sections on one slide.
- Avoid tears, cracks, or folds.
- Recommended section thickness for most tissue types is 10 µm.
- Tissues with higher fat content (e.g., breast tissue) may require sections closer to 20 µm.

Storage:

- Keep slides cold and transport on dry ice. **DO NOT** leave slides at room temperature.
- Slides need to be processed within 4 weeks so transfer them to Dalhousie Genomics Core as soon as possible after sectioning.
- Store slides individually (one slide per slide mailer or 50ml falcon tube) -80°C.

10X Genomics provides guidance on slide preparation for fresh frozen tissue which can be found [here](#).

FFPE Samples (Preferred)

RNA quality

- RNA quality of the tissue is assessed by calculating DV200 of RNA extracted from freshly collected tissue sections.
- Tissue with DV200 ≥30% are more likely to generate successful Visium Spatial Gene Expression libraries.
- If you'd like us to perform this QC step please supply two 20 µm FFPE sections prior to preparing your slides.

Section placement on slides:

- Place the tissue sections within the allowable area. Place one tissue section per slide unless the sections are small enough to be captured together.
- Ensure that tissue sections and associated paraffin **do not** overlap with other tissue sections and associated paraffin if placing multiple sections on one slide.
- Avoid tears, cracks, or folds.
- Recommended section thickness for most tissue types is 5 µm.

Storage:

- Maintain slides containing tissue sections in a low moisture environment such as a desiccator.
- Transport slides at room temperature.
- Slides need to be processed within 2 weeks so please ensure that you schedule sample submission before mounting tissue sections on slides.

10X Genomics provides guidance on slide preparation for FFPE tissue which can be found [here](#).

Fixed Frozen Samples**RNA quality**

- RNA quality assessment should be done prior to placing the tissue sections on slides.
- A DV200 score of ≥50% is recommended for Visium Gene Expression for Fixed Frozen samples.

Section placement on slides:

- Place the tissue sections within the allowable area. Place one tissue section per slide unless the sections are small enough to be captured together.
- Ensure that tissue sections and associated paraffin **do not** overlap with other tissue sections and associated paraffin if placing multiple sections on one slide.
- Avoid tears, cracks, or folds.
- Recommended section thickness for most tissue types is 10 µm.

Storage:

- Keep slides cold and transport on dry ice. **DO NOT** leave slides at room temperature.
- Slides need to be processed within 2 months.
- Store slides individually (one slide per slide mailer or 50-ml falcon tube) at -80°C.

10X Genomics provides guidance on slide preparation for Fixed Frozen samples which can be found [here](#).

Validated slides

| | Vendor | Part number |
|---|-------------------|--------------------|
| Epredia Shandon ColorFrost Plus Slides | Fisher Scientific | 6776214 |
| Fisherbrand Superfrost Plus Microscope Slides | Fisher Scientific | 12-550-15 |
| Poly-Prep Slides | Millipore Sigma | P0425 |
| VWR Superfrost Plus Micro Slides, Premium | VWR | 48311-703 |