


# NUCLEAR MAGNETIC RESONANCE RESEARCH RESOURCE

## Acquiring $^2\text{H}$ NMR Spectra

Mike Lumsden  
Coordinator, NMR-3  
August 16, 2021: version 2.00

### Before Beginning, Please Take Note of the Following...

- **Solvent:** Use non-deuterated solvents to avoid a  $^2\text{H}$  spectrum where the solvent signal towers over (and possibly obscures) your solute signals.
- **Locking & Shimming:**  $^2\text{H}$  NMR experiments must be acquired unlocked. In addition, a change in shimming procedure is required (instructions below).
- **Chemical Shift Referencing:** Since  $^2\text{H}$  NMR must be acquired unlocked, chemical shift referencing cannot be performed in the traditional manner. Options for referencing your  $^2\text{H}$  spectra include:
  - Use the Solvent  $^2\text{H}$  Signal: Use the natural abundance  $^2\text{H}$  signal or add a small amount of deuterated solvent to your protio solvent (~5%) for referencing purposes. The figure on the right shows the  $^2\text{H}$  spectrum obtained from tap water in 32 scans. This signal is from the isotopomer  $^1\text{H}^{16}\text{O}^2\text{H}$  (calculated concentration of 311 ppm or 16.3 mM).
  - Use the  $^1\text{H}$  Signal of TMS from a previously recorded proton spectrum: an au program exists making this approach trivial.
  - Use a separate NMR sample containing the deuterated version of your protio solvent.



### Instructions for $^2\text{H}$ NMR Experiments

- (1) **\*\*OPTIONAL\*\*** If you have a sample containing deuterated solvent for chemical shift referencing purposes, insert this sample first and lock. Once the instrument is locked, eject the sample.
- (2) Enter the command **“lockoff.ml”** to disable the lock channel. When complete (~ 10 s), an information window appears stating that the lock is off.
- (3) Insert your NMR sample in the magnet.

#### **\*\*\*CAUTION\*\*\***

Do not LOCK, SHIM or TUNE the PROBE while running  $^2\text{H}$  NMR Experiments!

- (4) Click NEW EXP in the button panel to create a new dataset for your  $^2\text{H}$  NMR experiment. Select **“1d\_2H”** from the experiment drop-down list.
- (5) Shim by typing the command **“lctshim.ml”**. This program will automatically perform a topshim using the most intense proton signal in your sample (usually from the solvent).
- (6) Change **“ns”** if desired (64 is the default).
- (7) Type **“gpro”** to ensure both the deuterium pulse width and power are up to date.
- (8) Click START in the button panel to acquire the  $^2\text{H}$  NMR data.
- (9) **\*\*OPTIONAL\*\*** Reference your spectrum using the  $^1\text{H}$  resonance frequency of TMS by typing the command **“xref.ml”**.
- (10) **\*\*OPTIONAL\*\*** If you have another sample for  $^2\text{H}$  NMR, return to step 3 above.
- (11) Eject your NMR sample from the magnet.
- (12) **\*\*IMPORTANT FINAL STEP\*\*** enter the command **“lockon.ml”** before finishing. This returns the lock channel back to the default configuration for the next user.