



The Vista Cervical Collar Inservice

Life Changing Spinal Innovation

Aspen Vista Collar Introduction

The Vista Cervical Collar from Aspen Medical Products is an adjustable collar. With its innovative height adjustment technology, the Vista is really six collars in one.



Life Changing Spinal Innovation

Aspen Vista Collar Introduction

By pulling out the sizing dial and turning clockwise, the Vista collar can be set at any of the six sizing options.




Upon completion of this program, the nurse will be able to:

- A. Describe the indications and usage of cervical collars
- B. Demonstrate the proper sizing and application of the cervical collar
- C. Demonstrate skin care protocols related to the cervical collar
- D. Demonstrate patient education related to the cervical collar

Aspen Vista Collar

The Vista Cervical collar consists of several pieces including the front panel, back panel, chin piece and collar pads. The Vista collar also includes a sizing dial and sizing indicator.



Aspen Vista Collar Packaging

The Vista collar is packaged at the lowest setting with the back panel wrapped around front panel. Prior to sizing, remove back panel and set aside.





Vista Collar Application



Life Changing Spinal Innovation

Vista Collar Application - Supine

There are 6 steps to fitting the Vista collar on a supine patient.

- A. Pre-form Collar**
- B. Position Back Panel**
- C. Position Front Panel**
- D. Adjust Height**
- E. Attach Back Panel**
- F. Check Fit**

Aspen Pre-Form Vista Collar

1. Pre-form the side panels by bending or rolling collar ends in toward center.
2. This will ensure that the collar conforms to patient anatomy and provides adequate motion restriction.



Aspen Position Back Panel

1. While holding the head in neutral alignment to maintain cervical spine alignment, Press down to compress bedding and slide back panel behind patient's neck.
2. The loop straps should be aligned between the patient's ear and trapezius muscle.



Aspen Position Front Panel

While pulling the side panels apart, position the collar under the patient's chin, in contact with the sternum, with the side panels up and over the trapezius muscles. This will prevent the collar from resting directly on the trapezius and causing discomfort.



Aspen Position Front Panel

1. At this point in the fitting process it is okay if the patient's chin is not touching the chin piece.
2. Ensure that the patient's chin is centered on the chin piece and aligned with the sternal notch.



Aspen. Adjust Height

While holding the front panel in place, grasp the dial between your thumb and forefinger. Note the indentations on either side of the side. Placing your thumb and forefinger will make it easier to pull the dial out slightly and release the lock.



Aspen. Adjust Height

Pull dial out and turn clockwise to raise chin piece.



Aspen. Adjust Height

Dial up until the foam pad on the chin piece begins to compress. It is important that the collar is properly sized to provide adequate motion restriction. Ensure there is not a gap between the patient's chin and the chin piece. When the chin piece is in proper position, release the dial and it will self-lock.



Aspen. Adjust Height

Ensure that the yellow dot is completely within the sizing indicator. After sizing, twist dial without pulling to verify it is locked in place.



Aspen Attach Back Panel

- A. While holding the front panel in place, attach the loop strap to the front panel on each side.
- B. Form the side panel to the patient's anatomy by pushing inward on the back edge while pulling the loop strap out. Tighten back panel by pulling loop strap up and reattach. Keep the back panel centered by tightening each side equally.



Aspen Attach Back Panel

- A. Repeat until all slack has been removed and a secure fit has been achieved.
- B. To check tightened collar, place forefinger inside collar, behind mandible on each side of patient. Simultaneously slide fingers forward towards chin. Collar should be tightened if you can easily slide your fingers all the way to the chin. This is to ensure that the collar is snug enough to provide adequate motion restriction.



Check to see that the back of the chin piece is not pressing inward on the patient's throat and the patient can breathe and swallow.



Aspen. Check Fit

If the chin piece is pressing on the patient's throat, loosen the back panel and lift the side panels up slightly, off of the trapezius muscles. This will tilt the front panel forward slightly, moving the chin piece forward and away from anterior aspect of the patient's throat.



Vista Collar Application - Upright

There are 6 steps to fitting the Vista collar on a supine patient.

- A. Pre-form Collar**
- B. Position Back Panel**
- C. Attach Front Panel**
- D. Adjust Height**
- E. Tighten Collar**
- F. Check Fit**

Aspen Pre-Form Vista[®] collar

1. Pre-form the side panels by bending or rolling collar ends in toward center.
2. This will ensure that the collar conforms to patient anatomy and provides adequate motion restriction.



Aspen Position Front Panel

While pulling the side panels apart, position the collar under the patient's chin, in contact with the sternum, with the side panels up and over the trapezius muscles. This will prevent the collar from resting directly on the trapezius and causing discomfort.



Aspen Position Front Panel

1. At this point in the fitting process it is okay if the patient's chin is not touching the chin piece.
2. Ensure that the patient's chin is centered on the chin piece and aligned with the sternal notch.



Aspen Attach Back Panel

1. Loosely attach back panel loop straps to white or gray areas to maintain front panel placement during sizing.
2. Do not attach loop strap to green area since this may interfere with the height adjustment of the front panel.



Aspen. Adjust Height

While holding the front panel in place, grasp the dial between your thumb and forefinger. Note the indentations on either side of the side. Placing your thumb and forefinger will make it easier to pull the dial out slightly and release the lock.



Aspen. Adjust Height

Pull dial out and turn clockwise to raise chin piece. When the chin piece is supporting the chin, release the height adjustment dial and it will self-lock. It is important that the collar is properly sized to provide adequate motion restriction. Ensure there is not a gap between the patient's chin and the chin piece. When the chin piece is in proper position, release the dial and it will self-lock.



Aspen Adjust Height

Ensure that the yellow dot is completely within the sizing indicator. After sizing, twist dial without pulling to verify it is locked in place.



Aspen[®] Tighten Collar

Form the side panel to the patient's anatomy by pushing inward on the back edge while pulling the loop strap out. Tighten back panel by pulling loop strap up and reattach. Keep the back panel centered by tightening each side equally.




Aspen Tighten Collar

- A. Repeat until all slack has been removed and a secure fit has been achieved.
- B. To check tightened collar, place forefinger inside collar, behind mandible on each side of patient. Simultaneously slide fingers forward towards chin. Collar should be tightened if you can easily slide your fingers all the way to the chin. This is to ensure that the collar is snug enough to provide adequate motion restriction.



Check to see that the back of the chin piece is not pressing inward on the patient's throat and the patient can breathe and swallow.



Aspen. Check Fit

If the chin piece is pressing on the patient's throat, loosen the back panel and lift the side panels up slightly, off of the trapezius muscles. This will tilt the front panel forward slightly, moving the chin piece forward and away from anterior aspect of the patient's throat.



Aspen Fitting Tip

Fitting Tip #1 – the patient's chin should be close to the front edge of the chin piece but the back of the chin piece should not be pressing on the patient's throat.



Aspen Fitting Tip

Fitting Tip #2 – all slack between the front and back panels must be removed to ensure a secure fit. If necessary, re-tighten both sides.



Fitting Tip #3 – patient's often change position, it may be necessary to readjust the height and the re-tighten.

For example, when a supine patient becomes ambulatory and gravity pulls down on the shoulders, it is often necessary to readjust the collar.

Fitting Tip #4 – if space is note beneath the chin piece, dial up the height. This will tighten the collar and may be all that is necessary for a proper fit. In some cases it may be necessary to tighten the loop straps evenly to achieve a proper fit.



Accessories – Back Panel Options

- A. The Vista collar comes with a regular back panel. For patients with larger anatomies, the Large Vista collar is available.
- B. For obese patients, the Extra Large Aspen back panel is available.





Replacement Pads



Life Changing Spinal Innovation

Aspen Changing Pads

The Vista pads can be removed and cleaned when they become soiled. Remove the pads from the collar by peeling away from hook dots.



Aspen. Cleaning the Collar

Wipe down the plastic shell of the collar using mild soap and water. Rinse and towel dry before reattaching pads.



Aspen. Cleaning the Pads

Hands wash pads with a mild soap and warm water.



Aspen. Cleaning the Pads

Rinse thoroughly to ensure no soap residue exists, this will prevent skin irritation.



Cleaning the Pads


- A. Allow pads to air dry. Do not use heat.
- B. Use replacement pads while first set is drying.

Aspen Vista® Cervical Collar Padding

The white side of the Vista padding is made of hypoallergenic cotton. This prevents skin irritation and effectively wicks moisture away from the skin. There are pads for the front panel, chin piece, and back panel.

Aspen. Replacing the Pads

To attach the chin piece pad to chin piece, fold pad so that white is touching white and press the gray side to the hook and loop dots.



Aspen. Replacing the Pads

To attach the front panel pad, fold pad so that white is touching white and press the gray side to the hook and loop dots. Center pad in the proper location on the plastic shell so that it extends past the edges of the collar. Make sure that the plastic collar frame does not come in contact with the patient's skin.



Aspen. Replacing the Pads

Place pad on back panel by pressing gray side to hook and loop dots. Thread loop strap through slots in pad.





Skin Care



Life Changing Spinal Innovation



Skin Care Protocol Recommendations

- A. Remove device at least once per shift.
- B. Make sure skin is clean and dry.
- C. Check for signs of skin breakdown.
- D. Palpate for skin changes within hairline.
- E. Change pads as necessary to ensure they are clean and dry.



Vista Collar Competency Measurement



Life Changing Spinal Innovation



Vista Collar competency tool



Life Changing Spinal Innovation

Vista Collar competency tool

Measurement 1 – Indications of cervical collars

Vista Collar competency tool

Measurement 1 – Usage of cervical collars

Vista Collar competency tool

Measurement 1 – Usage of cervical collars

Vista Collar competency tool

Measurement 2 – Correct measuring of Vista collar

Measurement 3 – Demonstrate proper application of collar

- A. Pre-form collar
- B. Position Back Panel
- C. Position Front Panel
- D. Sizing
- E. Attach and Tighten
- F. Check airway

Measurement 4 – Demonstrate Proper Skin Care Protocols

- A. Removal of Collar
- B. Skin Assessment and Cleaning
- C. Change Pads
- D. Reapply Collar

Measurement 5 – Patient Education

- A. Demonstrate how to remove pads
- B. Demonstrate how to wash pads
- C. Demonstrate how to re-apply pads to collar
- D. Provide patient with Vista® Cervical Collar Patient Handbook

Links

- A. Clinical Research
- B. Technical Reports
- C. Tools



Video Links



Life Changing Spinal Innovation

Test

Test



Slide for Large Images



Life Changing Spinal Innovation